



ecore™

Sustainable Solutions for Health, Productivity and the Environment.



# Enviro Care® Neutral Disinfectant



# Enviro Care® Neutral Disinfectant

## Where to Use

Recommended for use in hospitals, nursing homes, schools, colleges, commercial and industrial institutions, veterinary clinics, sports, fitness, and recreation applications. Disinfects, cleans and deodorizes floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, plastic surfaces and other hard non-porous surfaces.

### Benefits

- Hard Water Tolerance
- NPE Free
- Organic Soil Tolerance
- Neutral formula will not dull floor finishes
- Cleans, Disinfects, Deodorizes
- Pseudomonicidal, Fungicidal, Virucidal

### Dilution Ratios

Dilute: (1:64) 2 oz per gallon/60 ml per 3.8 litres of water.

### Directions for Use

**Application:** Remove heavy soil deposits from surface. Then thoroughly wet surface. The solution can be applied with a cloth, mop, sponge, or coarse spray, or soaking. Let solution remain on surface for a minimum of 10 minutes. Rinse or allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.

**For Human Coronavirus:** treated surfaces must remain wet for 1 minute. Wipe dry with a clean cloth, sponge or mop or allow to air dry.

**For Influenza Virus Type A:** treated surfaces must remain wet for 2 minutes. Wipe dry with a clean cloth, sponge or mop or allow to air dry.

Kills HIV-1 (AIDS virus) on pre-cleaned, environmental surfaces or objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human immunodeficiency virus Type I (HIV-1) (Associated with AIDS).

### Distributor Contact Info

Ecore International

Email: [salesupport@ecoreintl.com](mailto:salesupport@ecoreintl.com)

Phone: 800-322-1923

### Effective Against:

- HIV-1 (AIDS virus).
- Hepatitis B (HBV).
- Hepatitis C Virus (HCV).
- Staphylococcus aureus - Methicillin resistant strain (MRSA).
- Escherichia coli (antibiotic resistant).
- Aspergillus Niger, Candida albicans, Trichophyton mentagrophytes (a cause of ringworm).
- Pseudomonas aeruginosa (antibiotic resistant)
- Streptococcus faecalis (antibiotic resistant)
- Klebsiella pneumoniae (antibiotic resistant)

Properties	
Color	Blue
Fragrance	Mild
Viscosity	Water, thin
VOC's (% by weight) EPA Method 24	0.5%
pH (Concentrated)	8.2
Readily Biodegradable	Yes
Phosphate Free	Yes

Available in 1 gallon  
containers from Ecore.

1 gallon of Enviro Care makes 64 gallons  
of disinfectant.

For handling and precautionary information, reference  
the product label and safety data sheet.

For more information about  
RMC call 1.800.836.1633 or go  
to [www.rochestermidland.com](http://www.rochestermidland.com)

# Enviro Care<sup>®</sup>

# Neutral Disinfectant

## Efficacy Data

### Bactericidal

Enviro Care Neutral Disinfectant is bactericidal against many gram-negative pathogens according to the use of the dilution tests AOAC, modified in the presence of serum and hardness of organic water 5% 400 ppm (calculated in CaCO<sub>3</sub>).

Dilution 1:64 (2 oz per gallon / 60 mL for 3.8 liters of water).

Acinetobacter baumannii  
Acinetobacter calcoaceticus  
Bordetella bronchiseptica  
Chlamydia psittaci  
Enterobacter aerogenes  
Enterobacter cloacae  
Enterobacter cloacae -NDM-1  
Enterococcus faecalis (vancomycin resistant-VRE)  
Escherichia coli  
Escherichia coli-NDM-1  
Fusobacterium necrophorum  
Klebsiella pneumoniae  
Klebsiella pneumoniae-NDM-1  
Legionella pneumophila  
Listeria monocytogenes  
Pasteurella multocida  
Proteus mirabilis  
Proteus vulgaris  
Pseudomonas aeruginosa  
Salmonella choleraesuis  
Salmonella enterica  
Salmonella enteritidis  
Salmonella typhi  
Salmonella typhimurium  
Serratia marcescens  
Shigella flexneri  
Shigella sonnei

Community-associated methicillin-resistant Staphylococcus aureus - (CA-MRSA) (NRS384) (USA300)  
Community-associated methicillin-resistant Staphylococcus aureus - (CA-MRSA) (NRS123) (USA400)  
Staphylococcus aureus-résistant à la méthicilline (MRSA)  
Staphylococcus epidermidis  
Streptococcus pyogenes (Strep)  
Staphylococcus aureus-Vancomycin Intermediate Resistant (VISA)  
Streptococcus faecalis

### Virucidal

Enviro Care Neutral Disinfectant is virucidal against the following viruses after having been tested according to the virucide classification, according to the use of the dilution tests CaCO<sub>3</sub> modified in the presence of serum and hardness of the organic water 5% 400 ppm (calculated in CaCO<sub>3</sub>).

Dilution 1:64 (2 oz per gallon / 60 mL in 3.8 liters of water).

Adenovirus Type 7 - 8 oz. Per gallon  
Adenovirus Type 4  
Avian Influenza (H5N1)  
Avian polyomavirus  
Canine distemper

Coronavirus (SARS)- (cause of acute respiratory syndrome)  
Equine Herpes Virus Type 1  
Feline coronavirus  
Feline leukemia  
Feline panleukopenia virus  
Feline picornavirus  
Gastroenteritis transmissible (TGE)  
Hepatitis B Virus  
Hepatitis C virus  
Herpes Simplex Type I  
Herpes Simplex Type II  
HIV-1 (AIDS)  
Human coronavirus  
Infectious bovine rhinotracheitis  
Infectious bronchitis (Avian IBV)  
Influenza A / Hong Kong (H1N1)-Swine flu  
Newcastle disease  
Porcine parvovirus  
Pseudorabies (PRV)  
Rabies  
Respiratory syncytial virus (RSV)  
Rotavirus [WA]  
Rubella virus  
SARS Associated Coronavirus  
Vaccinia (Pox virus)

### Fungicidal

Enviro Care Neutral Disinfectant is effective against the following fungi according to the AOAC fungicide test modified in the presence of 5% iron (BIOSTRESS load) and water hardness of 400 ppm (CaCO<sub>3</sub> calculation) at 20°C.

Dilution 1:64 (2 oz per gallon / 60 ml in 3.8 liters of water).

Aspergillus niger  
Trichophyton mentagrophytes (athlete's foot)  
Candida albicans

### Sanitizing

Enviro Care Neutral Disinfectant is a sanitizer for hard nonporous non-food contact surfaces according to the method for Sanitizers - Non-Food Contact Surfaces, modified in the presence of 400 ppm hard water plus 5% organic serum in 60 seconds against:

Enterobacter aerogenes  
Staphylococcus aureus

TO SANITIZE NON-FOOD CONTACT HARD NON-POROUS SURFACES: Add 2 oz. of this product per gallon of water to sanitize hard nonporous surfaces. Treated surfaces must remain wet for 60 seconds. Wipe dry with a sponge, mop, or cloth or allow to air dry. Do not use in food preparation, food handling and food storage areas. This product is not to be used on glassware, dishware, cookware or utensils.